ORIGINAL ARTICLE

Evaluation of a conventional glass ionomer cement with new zinc formulation: effect of coating, aging and storage agents

Julius Zoergiebel · Nicoleta Ilie

Received: 27 October 2011 / Accepted: 9 April 2012 / Published online: 2 May 2012 © Springer-Verlag 2012

Abstract

Objective The study focused on a recently launched conventional glass ionomer cement (GIC) with a particular chemical formulation of both, filler and acrylic liquid, by analysing its mechanical behaviour in comparison to three conventional GICs. Furthermore, the effect of resin coating and storage conditions was evaluated.

Materials and methods Three commercially available GICs were chosen: Riva Self Cure (SDI), Fuji IX Fast (GC) and Fuji IX GP Extra/Equia (GC). Additionally a newly developed zinc-containing GIC—ChemFil Rock (Dentsply)—was tested. Mechanical properties were determined at macro- [flexural strength (FS) and modulus of elasticity (E_{flexural})] and microscale [Vickers hardness (VH) and indentation modulus (E)] after storing coated and uncoated specimens in artificial saliva and distilled water for 7 and 30 days.

Results ChemFil Rock revealed the highest FS, but the lowest VH and *E*. The micro-mechanical properties of the analysed GICs did neither benefit from the new zinc formulation nor from resin coating. A resin coating is nevertheless a valuable support for GIC fillings, since it offers the absence of visible surface defects like crazing and voids, and thus, it led to significant improvements in flexural strength. This statement is also valid for ChemFil Rock, contrary to manufacture recommendation. The impact of storage agent and storage duration on the measured properties was low. *Conclusions* The new development (ChemFil Rock) might represent a promising approach regarding longevity of GIC

J. Zoergiebel \cdot N. Ilie (\boxtimes)

Department of Restorative Dentistry, Dental School of the Ludwig-Maximilians-University, Goethestr, 70, 80336, Munich, Germany e-mail: nilie@dent.med.uni-muenchen.de fillings in molar regions, due to the high flexural strength and the absence of visible surface defects like crazing and voids. *Clinical relevance* All GICs should receive surface protection in order to perform their maximum in stability.

Keywords Glass ionomer cements · Coating · Storage · Micro-mechanical properties

Introduction

In the field of dentistry, glass ionomer cements (GICs) are a common and useful choice for restorative therapy concerning fillings which are not situated in high-stress sites. However compared to permanent filling materials like resin-based composites, GICs show several advantages, such as the ability to adhere to moist enamel and dentin without necessitating an intermediate agent and anti-cariogenic properties such as the long-term fluoride release. Other clinical advantages like biocompatibility and low coefficient of thermal expansion support their valuable position in the daily dental practice [1-5].

These positive properties are unfortunately dwarfed by a poor surface polish, a high porosity and rather weak mechanical properties, such as brittleness, surface wear or fracture toughness [5, 6]. Therefore it is doubtful that GIC represents a capable counterpart of amalgam or resin-based composites in high-stress sites.

The first 10 min of the hardening process of GICs is characterised by a slow release of calcium ions within the matrix, followed by aluminium ions [7]. During this time period, the material is very much frail to dehydration [5], meaning in clinical conditions that a gain or loss of liquid in this phase can tremendously affect the final properties of the restoration. One advance in strengthening GICs was the application of surface protection in order to preserve the balance of the amount of liquid [8–11]. Reviewing the literature, only a few studies are dealing with the effect of resin coating on conventional GICs. Earl et al. [8] found that immediate covering of the immature GIC surface with light-activated bonding resin is the most effective method of limiting water movement across the surface. This supports the balance between water uptake and loss, leading to better clinical results. In a study by Ribeiro et al. [10], in analysing the effectiveness of surface protection for resin-modified GICs, it was stated that all tested materials required surface protection. Furthermore they found that Heliobond light-activated bonding resin displayed superior results than nail varnishes and surface coatings suggested by the manufacturer.

The manufacturer of a recently launched GIC (ChemFil Rock, Dentsply) followed a different approach to enhance material's stability, claiming that surface protection in the form of resin coating is irrelevant for product's performance. An enhanced setting reaction in the new GIC is supposed, due to the zinc content as part of its glass particles, leading thus to higher strength, by similar working time and application comfort as regular GICs [12].

Our present study aimed therefore primarily to verify if the zinc-containing material performed comparable to traditional alumina–silicate GICs. The influence of resin coating, aging and aging agent, such as artificial saliva and distilled water, are evaluated. The null hypotheses tested were that: (a) the new GIC would perform similar to the traditional GIC in terms of macro- (flexural strength and modulus of elasticity in flexural test) and micro- (Vickers hardness and modulus of elasticity) mechanical properties, (b) applying resin coating on the GICs surface would not influence the above-mentioned properties, (c) aging (7 and 30 days) would not affect the measured mechanical properties and (d) storing the GICs in artificial saliva or distilled water would result in similar mechanical properties.

Materials and methods

Three commercially available conventional restorative GICs— Riva Self Cure (SC), Fuji IX GP Fast and Fuji IX GP Extra (Equia)—and a new conventional GIC, ChemFil Rock, were selected (Table 1). Additionally, the corresponding light-cured resin coating for each material was chosen (Table 1). It should be noticed that for ChemFil Rock, no surface protection is indicated by the manufacturer. We still choose to apply an experimental resin coating supplied by the same manufacturer in order to confront every material with the same conditions. The mechanical properties were determined at macro- [flexural strength (FS) and modulus of elasticity ($E_{flexural}$)] and microscale [Vickers hardness (VH) and indentation modulus (E)] for coated and uncoated specimens after storing the samples in artificial saliva or distilled water for 7 and 30 days, respectively. In order to evaluate the size and shape of the glass particles, a scanning electron microscope (SEM) was used.

Macro-mechanical characteristics

Using a three-point bending test (in analogy to ISO 4049:2009) on bar-shaped specimens ($16 \times 2 \times 2$ mm), FS and E_{flexural} (n=20) were evaluated. Four materials (Table 1), two coating conditions (with and without coating), two aging times (7 and 30 days) and two aging agents (distilled water and artificial saliva) resulted in 32 groups, summarising 640 samples.

The encapsulated GICs were mixed by rotating in a RotoMix (3M-ESPE, Seefeld, Germany) apparatus and fabricated at room temperature, according to the manufacturer instructions. The specimens were prepared in a stainless steel mould, allowing them to set for 20 min at room temperature in the mould. In order to achieve plane and possibly voids free specimens, a transparent foil (US-120 KE; Frasaco, Tettnang, Germany) was placed on the surface of the unset material, by pressing it with an object slide. The foils were removed about 2.5 min after mixing. Specimens, receiving resin coating, were prepared in the same way, except for the application of the coating. The latter was applied 3 min after activation of the GIC capsules. Thereby the transparent foils were removed and the coating product was applied. Another foil was gently pressed on the coated surface of the specimen followed by light curing the coating for 20 s (Mini L.E.D, SATELEC SED-R, France) by three overlapping irradiations. After 20 min the specimens were taken out of the moulds and conditioned either in artificial saliva or in distilled water at 37°C, and stored for 7 or 30 days.

Previously to loading into a universal testing machine (MCE 2000ST; Quicktest Prüfpartner GmbH, Langenfeld, Germany), each specimen was gently grounded with 1,200grit silicon carbide (SiC) paper (LECO, St. Joseph, MI, USA) and the exact dimension was recorded. The load was applied at a constant crosshead speed of 0.5 mm/min until fracture. The distance between the supporting points was 12 mm. During loading, the specimens were immersed in distilled water at room temperature. The coated side of the specimen was facing towards the tensile zone. E_{flexural} was calculated from the linear part of the force–deflection diagram.

Micro-mechanical characteristics

Fragments of the three-point bending test were used to determine VH and E by means of an automatic micro-hardness indenter (Fischerscope H100C, Fischer, Germany). Ten randomly selected specimens of each group were wet grounded with 2,500 and 4,000-grit SiC paper (FEPA). The coating was
 Table 1
 Materials, manufacturer and chemical composition of glass ionomer cements (all encapsulated) and coating materials

	Manufacturer	Composition
Glass ionomer cement		
Riva Self Cure;	SDI Limited,	Fluoro-aluminosilicate glass
Lot: B1004281	Victoria, Australia	Polyacrylic acid + tartaric acid, polyacrylic acid
ChemFil Rock; Lot: 1005004004	Dentsply, Konstanz, Germany	Calcium-aluminium-zinc-fluoro-phosphor-silicate glass, polycarboxylic acid, iron oxide pigments, titanium dioxide pigments, tartaric acid, water
GC Fuji IX GP Fast; Lot: 1005211	GC Europe N.V., Leuven, Belgium	Alumino-fluoro-silicate glass, polyacrylic acid, distilled water, polybasic carboxylic acid
GC Fuji IX GP Extra (Equia); Lot: 1005281	GC Europe N.V., Leuven, Belgium	Polyacrylic acid, aluminosilicate glass, distilled water
Coating		
Riva Coat; Lot: 091103	SDI Limited, Victoria, Australia	Acrylic monomer
Seal&Protect TF; Lot: MTO-3-27-1	Dentsply, Konstanz, Germany	Di- and trimethacrylate, acetone, dipentaerythritol penta acrylate monophosphate
GC Fuji Coat LC; Lot: 1005061	GC Europe N.V., Leuven, Belgium	Methylmethacrylate, multifunctional methacrylate, camphorquinone
GC G-Coat Plus; Lot:1004091	GC Europe N.V., Leuven, Belgium	Methyl methacrylate, colloidal silica, camphorquinone

completely removed by this procedure. Six indentations were accomplished on each sample's coated sight, with a total of 60 measurements in each group. The measurements were carried out force controlled: The test load increased and decreased with constant speed between 0.4 and 500 mN. The load and the penetration depth of the indenter were continuously measured during the load–unload hysteresis.

The universal hardness is defined as the test force divided by the apparent area of the indentation under the applied test force. From a multiplicity of measurements stored in a database supplied by the manufacturer, a conversion factor between universal hardness and Vickers hardness was calculated and implemented into the software, so that the measurement results were indicated in the more familiar Vickers hardness. The indentation modulus was calculated from the slope of the tangent of indentation depth curve at maximum force (DIN-50359-1 (1997) Testing of metallic materials—universal hardness test—part 1: test method).

The samples used for the SEM analysis (Zeiss Supra 55 VP, Oberkochen, Germany) were stored for 7 days in distilled water and received no sputtering. The images were taken by using a backscatter signal (RBSD).

Statistical analysis

The Kolmogorov–Smirnov test was applied to verify if the data were normally distributed. Results were compared using one and multiple-way ANOVA and Tukey's HSD post hoc test (α =0.05). An independent *t* test additionally analysed the differences in mechanical properties as function of coating, aging duration and aging agent (SPSS Inc.;

Chicago, IL, USA, Version 19.0). An additional Weibull analysis was performed for the flexural strength data. A multivariate analysis (general linear model with partial etasquared statistics) assessed the effect's strength of the parameters GIC, coating, storage agent and storage duration on the considered properties. The partial eta-squared statistic reports the practical significance of each term, based upon the ratio of the variation accounted for by the effect. Larger values of partial eta squared indicate a greater amount of variation accounted for by the model effect, to a maximum of 1.

Results

Post hoc multiple pairwise comparisons with Tukey's HSD test (p < 0.05) revealed the highest values in FS for ChemFil Rock under all measured conditions (Tables 2 and 3). The influence of coating on the flexural strength was proved to be significant (p < 0.05) for almost all tested conditions (exception: Riva Self Cure, saliva, 1 month and Riva Self Cure water, 1 week). Furthermore it could be verified that aging has no significant effect on FS (exception: Riva Self Cure with coating, stored in water) in both of the analysed storage agents. As for the Weibull parameter *m*, the reliability of the material, no clear dependency from material, coating condition, aging agent or duration can be emphasized. The highest value (*m*) was generated in ChemFil Rock, coated and stored for 30 days in artificial saliva.

Regarding the flexural modulus, E_{flexural} , the difference among materials was lower. Similar is valid also for the influence of the storage agent. A significant higher flexural

Material	Storage	Coating					No coating					Saliva	Water
		Saliva	ш	Water	ш	p1	Saliva	ш	Water	ш	p2	p3	p4
Riva Self Cure	1W	19.3a (4.9)	4.32	14.3A (8.3)	1.73	0.033	13.8a (4.4)	4.12	13.3A (3.7)	4.29	0.73	0.001	0.641
	1M	17.6α (9.9)	1.52	22.8B (7.8)	2.47	0.072	15.5a (5.8)	3.22	12.8A (3.4)	4.64	0.09	0.414	0
ChemFil Rock	1W	38.7ß (12.5)	3.77	39.1C (7.7)	5.66	0.908	28.2b (11.1)	2.74	31.7B (9.6)	3.15	0.308	0.008	0.011
	1M	41.5β (3.8)	13.1	39.4C (6.6)	7.08	0.22	30.9b (8.2)	4.19	30.4B (10.5)	3.23	0.867	0	0.003
Fuji IX GP Fast	1W	19.4α (6.2)	3.74	16.7AB (5.9)	4.39	0.17	12.7a (8.2)	1.97	9.1A (1.6)	6.88	0.079	0.007	0
	1M	23.0α (4.7)	5.76	15.5A (4.3)	4.59	0	9.7a (2.6)	4.17	9.2A (1.9)	5.53	0.539	0	0
Equia	1W	20.4α (3.4)	6.81	20.1AB (4.9)	4.63	0.815	12.0a (7.4)	2.21	9.7A (4.1)	3.04	0.221	0	0
	1M	22.4α (4.7)	5.83	18.5AB (5.8)	3.63	0.027	9.2a (5.0)	2.51	12.6A (6.9)	2.13	0.1	0	0.007

m is indicated

in saliva (p3) and in water (p4)). The Weibull parameter

no coating for samples stored

vs.

samples with coating (p1) and without coating (p2)) and coating (difference coating

modulus (p<0.05) when samples were stored in artificial saliva compared to samples stored in water was found only in few groups (both, coated and uncoated samples of Riva Self Cure, 1 week and Fuji IX GP Fast, 1 month of storage). The positive influence of coating on the modulus of elasticity was more frequent evident in samples stored in water (Fuji IX GP Fast and Equia all storage durations) than in saliva (Riva Self Cure and Equia, both materials 1 week of storage).

The micro-mechanical properties, in contrast, varied consistently among the tested materials [Table 4 (a and b)], achieving the highest VH in Fuji IX Fast with coating, stored for 4 weeks in water. Generally this GIC reached the significant highest VH and E values under all measured conditions, whereas ChemFil Rock the significant lowest.

A significant increase in both measured micromechanical properties, VH and *E*, with storage duration was measured only for Riva Self Cure stored in saliva and uncoated. The influence of coating and storage agent on the micro-mechanical properties was in most cases not significant. Considering the multivariate analysis (Table 5), the flexural strength was proved to be stronger influenced by the material itself (η^2 =0.915) than by the coating condition (η^2 =0.740), whereas the effect on the flexural modulus was consistently lower (lower η^2 values, Table 5), with the coating condition (η^2 =0.236) exerting a stronger influence than the storage duration (η^2 =0.190).

The SEM images display small glass particles size in ChemFil Rock (Fig. 1). The diameter of the particles differs consistently from the larger diameter of the other three measured materials, demonstrating a greater particle surface area in ChemFil Rock. The glass particles in ChemFil Rock appeared rather homogenous in size and shape (SEM images, Fig. 1). More erratically glass particles were found in Riva SC, Fuji IX Fast and Equia.

Discussion

The study evaluated whether a new developed GIC (Chem-Fil Rock, Dentsply) can provide higher mechanical properties than conventional GIC containing alumina silicate glasses. According to the manufacturer, this new material should present an enhanced durability due to a zinccontaining glass with an accelerated ion release pattern, when compared to conventional GICs. Besides changes in the chemical composition of the glass particle, also a novel acrylic acid copolymer with a high molecular weight was incorporated in the material [12].

Generally, the setting process in a GIC is characterised by an interaction between a polyacid liquid and a glass powder in form of an acid–base reaction [13]. Upon an initial buildup of calcium polyalkenoate, the formation of aluminium

Material	Storage	Coating	Coating			No Coating			
		Saliva	Water	p1	Saliva	Water	p2	p3	p4
Riva Self Cure	1W	6.3αβ (1.3)	4.3A (1.3)	0	5.4a (0.8)	4.4A (0.9)	0.001	0.016	0.72
	1M	5.4α (1.7)	6.4B (1.5)	0.062	5.5a (2.5)	5.8BC (0.8)	0.538	0.933	0.167
ChemFil Rock	1W	5.9αβ (2.1)	5.9B (1.1)	0.919	5.3a (1.0)	5.3AB (1.0)	0.979	0.286	0.141
	1M	6.4αβ (1.1)	6.1B (1.0)	0.367	6.0a (1.2)	6.6C (1.0)	0.072	0.203	0.133
Fuji IX GP Fast	1W	6.6αβ (2.0)	7.0B (2.3)	0.548	5.9a (1.5)	5.8BC (0.8)	0.878	0.238	0.045
	1M	7.2β (2.0)	6.0B (1.1)	0.02	6.4a (1.0)	5.0AB (1.1)	0	0.134	0.011
Equia	1W	6.3αβ (1.0)	6.1B (1.0)	0.403	5.3a (1.4)	5.0AB (1.0)	0.489	0.01	0.003
-	1M	6.5αβ (1.1)	6.6B (1.0)	0.74	5.9a (1.1)	5.4B (1.2)	0.176	0.145	0.004

Table 3 Macro-mechanical properties of modulus of elasticity in flexural test (E_{flexural} , in gigapascals) are detailed in mean values and standard deviations (in parentheses)

Symbols/letters indicate statistically homogeneous subgroups within a column (Tukey's HSD test, α =0.05). A *t* test analysed differences as function of storage (difference saliva vs. water storage for samples with coating (p1) and without coating (p2)) and coating (difference coating vs. no coating for samples stored in saliva (p3) and in water (p4)). The Weibull parameter *m* is indicated

polyalkenoate occurs. This reaction is characterised by a stepwise, rather long-lasting setting, where changes in mechanical properties occur, mainly characterised by a rise in strength within the first 24 h. A continuing altering of strength can be observed over several weeks and months [14–16].

At first, the precipitation of the cement is continuing until most of the ions are in insoluble form. This period can be observed 3 to 6 min after the mixing process [17, 18].

Considering clinical conditions, both, water contamination and dehydration in this phase can be compromising for the cement applied in a cavity [5]. Crack propagation, frequently observed in GIC fillings, is a typical mechanical failure resulting from desiccation [19]. Gemalmaz et al. [20] observed that due to early moisture contamination, the mechanical properties of GICs decreased and their surfaces became more susceptible for erosion and abrasion. Naasan and Watson [5] emphasizes that water contamination must

 Table 4
 Micro-mechanical property of (a) Vickers hardness and (b) indentation modulus is detailed in mean values and standard deviations (in parentheses)

Material	Storage	Coating			No Coating			Saliva	Water
		Saliva	Water	p1	Saliva	Water	p2	p3	p4
Riva Self Cure	1W	105.6βγ (7.4)	116.1BC (7.4)	0.034	85.3bc (6.1)	101.8CD (33.8)	0.289	0	0.354
	1M	114.5βγδ (11.4)	115.5BC (11.3)	0.886	111.5de (23.6)	92.5BC (9.6)	0.112	0.789	0.004
ChemFil Rock	1W	59.1a (5.3)	67.6A (5.9)	0.024	56.4a (3.0)	60.8A (2.1)	0.016	0.313	0.036
	1M	68.4a (5.1)	66.8A (1.9)	0.498	64.8ab (2.6)	67.1AB (2.2)	0.127	0.161	0.814
Fuji IX GP Fast	1W	119.3γδ (3.0)	119.3BC (11.3)	0.998	113.0de (6.7)	118.5CD (13.9)	0.406	0.073	0.921
	1M	124.68 (11.9)	125.9C (19.0)	0.896	115.0e (14.6)	122.1D (10.0)	0.349	0.24	0.683
Equia	1W	103.0β (9.2)	99.3B (14.3)	0.603	90.7cd (16.5)	97.6CD (7.5)	0.383	0.149	0.802
	1M	110.9βγδ (5.2)	109.7BC (13.8)	0.846	97.9cde (10.4)	103.3CD (11.8)	0.423	0.027	0.407
Riva Self Cure	1W	24.1δε (1.1)	24.3CD (1.3)	0.746	19.6b (2.5)	22.0CD (2.2)	0.119	0.003	0.061
	1M	22.5γδ (1.6)	23.1C (0.8)	0.437	23.1c (1.7)	20.5BC (1.5)	0.018	0.571	0.006
ChemFil Rock	1W	17.4α (0.6)	18.9AB (0.9)	0.008	17.0a (0.6)	17.8A (0.7)	0.072	0.31	0.039
	1M	16.5a (1.3)	17.1A (0.8)	0.34	17.1a (0.4)	17.2A (0.45)	0.907	0.279	0.856
GC Fuji IX GP Fast	1W	24.3δε (0.7)	25.6D (1.7)	0.129	23.9c (0.6)	23.7DE (1.2)	0.642	0.291	0.044
	1M	25.3ε (1.1)	25.7D (1.7)	0.657	24.5c (1.7)	24.4E (1.2)	0.87	0.371	0.159
Equia	1W	19.9β (1.1)	18.7AB (1.2)	0.115	18.1ab (0.8)	18.7AB (0.8)	0.189	0.01	0.977
	1M	20.9βγ (1.0)	20.0B (1.2)	0.218	19.3ab (0.8)	19.4AB (0.6)	0.885	0.018	0.276
	Material Riva Self Cure ChemFil Rock Fuji IX GP Fast Equia Riva Self Cure ChemFil Rock GC Fuji IX GP Fast Equia	Material Storage Riva Self Cure IW ChemFil Rock IW Fuji IX GP Fast IW Riva Self Cure IW Riva Self Cure IW ChemFil Rock IW AM ChemFil Rock IW IM Sub IW IM IM IM IM IM IM IM IM IM IM	Material Storage Coating Riva Self Cure 1W 105.6βγ (7.4) 1M 114.5βγδ (11.4) ChemFil Rock 1W 59.1α (5.3) 1M 68.4α (5.1) Fuji IX GP Fast 1W 119.3γδ (3.0) 1M 124.6δ (11.9) Equia 1W 103.0β (9.2) Riva Self Cure 1W 103.0β (9.2) Riva Self Cure 1W 24.1δε (1.1) 1M 22.5γδ (1.6) 1 ChemFil Rock 1W 17.4α (0.6) IM 16.5α (1.3) 16.5α (1.3) GC Fuji IX 1W 24.3δε (0.7) GP Fast 1M 25.3ε (1.1) Equia 1W 19.9β (1.1)	Material Storage Coating Riva Self Cure IW 105.6βγ (7.4) 116.1BC (7.4) IM 114.5βγδ (11.4) 115.5BC (11.3) ChemFil Rock IW 59.1α (5.3) 67.6A (5.9) IM 68.4α (5.1) 66.8A (1.9) Fuji IX GP Fast IW 119.3γδ (3.0) 119.3BC (11.3) Fuji IX GP Fast IW 119.3γδ (3.0) 125.9C (19.0) Equia IW 103.0β (9.2) 99.3B (14.3) IM 109βγδ (5.2) 109.7BC (13.8) Riva Self Cure IW 24.1δε (1.1) 24.3CD (1.3) IM 22.5γδ (1.6) 23.1C (0.8) ChemFil Rock IW 17.4α (0.6) 18.9AB (0.9) IM 16.5α (1.3) 17.1A (0.8) GC Fuji IX GP Fast IW 24.3δε (0.7) 25.6D (1.7) GP Fast IM 25.3ε (1.1) 25.7D (1.7) Equia IW 19.9β (1.1) 18.7AB (1.2)	Material Storage Coating Riva Self Cure 1W 105.6βγ (7.4) 116.1BC (7.4) 0.034 1M 114.5βγδ (11.4) 115.5BC (11.3) 0.886 ChemFil Rock 1W 59.1α (5.3) 67.6A (5.9) 0.024 1M 68.4α (5.1) 66.8A (1.9) 0.498 Fuji IX GP Fast 1W 119.3γδ (3.0) 119.3BC (11.3) 0.998 Equia 1W 103.0β (9.2) 99.3B (14.3) 0.603 IM 104.0β (5.2) 109.7BC (13.8) 0.846 Riva Self Cure 1W 24.3δε (1.1) 24.3CD (1.3) 0.746 IM 16.5α (1.3) 17.1A (0.8)	MaterialStorageCoatingNo CoatingRiva Self CureIW105.6βγ (7.4)I16.1BC (7.4)p1SalivaRiva Self CureIW105.6βγ (7.4)116.1BC (7.4)0.03485.3bc (6.1)IM114.5βγδ (11.4)115.5BC (11.3)0.886111.5de (23.6)ChemFil RockIW59.1α (5.3)67.6A (5.9)0.02456.4a (3.0)IM68.4α (5.1)66.8A (1.9)0.49864.8ab (2.6)Fuji IX GP FastIW119.3γδ (3.0)119.3BC (11.3)0.998113.0de (6.7)IM124.6δ (11.9)125.9C (19.0)0.896115.0e (14.6)EquiaIW103.0β (9.2)99.3B (14.3)0.60390.7cd (16.5)IM110.9βγδ (5.2)109.7BC (13.8)0.84697.9cde (10.4)Riva Self CureIW24.1δε (1.1)24.3CD (1.3)0.74619.6b (2.5)IM22.5γδ (1.6)23.1C (0.8)0.43723.1c (1.7)ChemFil RockIW17.4α (0.6)18.9AB (0.9)0.00817.0a (0.6)IM16.5α (1.3)17.1A (0.8)0.3417.1a (0.4)GC Fuji IX GP FastIM25.3ε (1.1)25.7D (1.7)0.65724.5c (1.7)EquiaIW19.9β (1.1)18.7AB (1.2)0.11518.1ab (0.8)IM20.9βγ (1.0)20.0B (1.2)0.21819.3ab (0.8)	MaterialStorageCoatingNo CoatingRiva Self CureIWIo5.6βγ (7.4)I16.1BC (7.4)p1SalivaWaterRiva Self CureIW105.6βγ (7.4)I16.1BC (7.4)0.03485.3bc (6.1)101.8CD (33.8)IM114.5βγδ (11.4)115.5BC (11.3)0.886111.5de (23.6)92.5BC (9.6)ChemFil RockIW59.1a (5.3)67.6A (5.9)0.02456.4a (3.0)60.8A (2.1)Fuji IX GP FastIW119.3γδ (3.0)119.3BC (11.3)0.998113.0de (6.7)118.5CD (13.9)EquiaIW119.3γβ (3.0)119.3BC (11.3)0.908115.0e (14.6)122.1D (10.0)EquiaIW103.0β (9.2)99.3B (14.3)0.60390.7cd (16.5)97.6CD (7.5)Riva Self CureIW24.1δε (1.1)24.3CD (1.3)0.46497.9cde (10.4)103.3CD (11.8)Riva Self CureIW24.1δε (1.1)24.3CD (1.3)0.74619.6b (2.5)22.0CD (2.2)IM22.5γδ (1.6)23.1C (0.8)0.34417.1a (0.4)17.2A (0.45)GC Fuji IX GP FastIW24.3δε (0.7)25.6D (1.7)0.1223.9c (0.6)23.7DE (1.2)EquiaIW19.9β (1.1)18.7AB (1.2)0.11518.1ab (0.8)18.7AB (0.8)IM20.9βγ (1.0)20.0B (1.2)0.21819.3ab (0.8)19.4AB (0.6)	MaterialStorageCoatingNo CoatingRiva Self CureIW105.6βγ (7.4)116.1BC (7.4)0.03485.3bc (6.1)101.8CD (3.3.8)0.289IM114.5βγδ (11.4)115.5BC (11.3)0.886111.5de (23.6)92.5BC (9.6)0.112ChemFil RockIW59.1α (5.3)67.6A (5.9)0.02456.4a (3.0)60.8A (2.1)0.161Fuji IX GP FastIW119.3γδ (3.0)119.3BC (11.3)0.998113.0de (6.7)118.5CD (13.9)0.404Fuji IX GP FastIW119.3γδ (3.0)119.3BC (11.3)0.908115.0e (14.6)122.1D (10.0)0.349Fuji IX GP FastIW103.0β (9.2)99.3B (14.3)0.60390.7cd (16.5)97.6CD (7.5)0.383Riva Self CureIW103.0β (9.2)109.7BC (13.8)0.48497.9cde (10.4)103.3CD (11.8)0.424Riva Self CureIW24.1δε (1.1)24.3CD (1.3)0.7419.6b (2.5)22.0CD (2.2)0.119Riva Self CureIW24.1δε (1.1)24.3CD (1.3)0.7419.6b (2.5)22.0CD (2.2)0.119GC Fuji IX GP FastIW17.4α (0.6)18.9AB (0.9)0.0817.1a (0.4)17.2A (0.45)0.907GC Fuji IX GP FastIW25.3ε (1.1)25.7D (1.7)0.65724.5c (1.7)24.4E (1.2)0.87FujiIM25.3ε (1.1)25.7D (1.7)0.65724.5c (1.7)24.4E (1.2)0.87FujiIM20.9βγ (1.0)20.0B (1.2)0.11518.1ab (0.8) <t< td=""><td>Material Storage Coating No Coating Material No Coating Material Saliva Material No Coating Material Saliva Material Material Material Saliva Material Material Saliva Material Material Saliva Saliva</td></t<>	Material Storage Coating No Coating Material No Coating Material Saliva Material No Coating Material Saliva Material Material Material Saliva Material Material Saliva Material Material Saliva Saliva

Symbols/letters indicate statistically homogeneous subgroups within a column (Tukey's HSD test, α =0.05). A *t* test analysed differences as function of storage (difference saliva vs. water storage for samples with coating (p1) and without coating (p2)) and coating (difference coating vs. no coating for samples stored in saliva (p3) and in water (p4))

Table 5 Influence of material, storage agent, storage duration and coating on flexural strength (FS), modulus of elasticity in flexural test (E_{flexural}), Vickers hardness (VH) and indentation modulus (*E*)

Parameter	Eflexural	FS	VH	Ε
Material	NS	0.915	0.775	0.911
Storage agent	NS	NS	NS	NS
Storage duration	0.190	NS	NS	NS
Coating	0.236	0.740	NS	NS
Material Storage agent Storage duration Coating	NS NS 0.190 0.236	0.915 NS NS 0.740	0.775 NS NS NS	0.9 NS NS NS

The higher the partial eta-squared values, the higher is the influence of the selected variables on the measured properties [general linear model (η^2)]

NS nonsignificant

be prevented during the delicate phase of the setting, which is referred by the author to endure 1 day up to 2 weeks. The same author [5] also claims that as time progresses, surface coatings fade away on the strength of mastication wear. Within this process, the resistance of the cement increases towards variations in water balance due to post hardening [5].

Regarding the results of our study, it becomes obvious that once the materials were protected from the influence of aqueous solutions by light-cured resin, FS could be improved. Hereby, it has to be pointed out that all samples prepared in this study were protected by a transparent foil in the initial phase of setting, which in this case was about 150 s. Thus, an initial surface protection was actually supplied for all specimens, coated and uncoated, which cannot be directly transcribed to real clinical conditions. This step was however necessary to allow correct sample preparation for the three-point flexural strength test, which is seen as the most meaningful test to evaluate the mechanical behaviour of GICs [21, 22]. As described in former investigations [21], FS shows a high sensitivity towards surface irregularities, resulting in erosions caused by water. Crisp and Wilson [7] presumed that the high affinity of water to GICs is caused by the ion-depleted siliceous phase, whose behaviour is analogical to silica gel, taking up water from the surroundings. Xie et al. [23] conclude that a less dense surface, or rather larger and higher amounts of voids result in worse mechanical properties. The assumption can be made that FS usually rises when surface protection is applied. But the strength of GIC is deciding influenced by its glass composition as well [24]. Especially a high content of fluoride was proved to induce a higher compressive and flexural strength [24]. Similar is valid also for GIC containing high amounts of zinc embodied in the glass powder, since an enhanced network connectivity will occur, thus raising the ability of the material to form a cement with acrylic acid. The setting time was also shown to decrease, making the resulting GIC more resistant against hydrolysis and, finally, inducing in the material a higher strength [25]. The high FS measured for ChemFil Rock compared to the other materials in this study supports this thesis. Further examinations upon the cross-link density in GIC by using,

Fig. 1 Scanning electron microscopy images (backscattered modus on not sputtered samples). a ChemFil Rock, b Riva Self Cure, c Fuji IX GP Fast and d Equia



among others, MAS-NMR spectroscopy figured out that zinc oxide is primarily a network modifier than a network former [26]. It was also evidenced that an increasing amount of zinc is directly connected to an enhanced reactivity [27]. As a network modifier, zinc oxide contributes to the Si–O– Si bond disruption in the glass and thus increases the vulnerability of the glasses to acid attack. According to the manufacturer [12] and confirmed by our data, the formation of zinc–polycarboxylate complexes during the setting of the novel GIC enhances the strengths more than other complexes consisting of bivalent ions like calcium or strontium.

Besides the chemical composition of glasses, also the polyacrylate acids in GIC can influence the setting reaction and the resulting properties of GIC [27]. In ChemFil Rock, a new acrylic acid copolymer with a high molecular weight was incorporated as well. The material also has included in its chemical composition itaconic acid as a co-monomer, which when incorporated in a conventional commercial GIC was shown to improve the biaxial flexural strength and diametral tensile strength compared to compositions without this copolymer [28]. The itaconic acid as a co-monomer to the high molecular polycarboxylic acid contained in ChemFil Rock is supposed to reduce the interaction between the high molecular polyacids in form of hydrogen bridges, and thus to delay the building of a gel phase which would worsen the storage stability [12]. The working properties of composition with or without itaconic acid as a co-monomer were proved to be comparable and acceptable for water-based cements [28].

Another approach to explain the high flexural strength of ChemFil Rock is suggested by Prentice et al. [29] who found that improved strength of GIC is related to a decrease in mean particle size and thus an increase in glass surface area. This thesis might correspond to our results within the SEM imaging. Moreover, the visual inspection of the GIC samples previously to loading in the three-point bending showed crazing and surface imperfections for the materials Riva Self Cure, Fuji IX Fast and Equia, but none for ChemFil Rock. There is thus evidence that macro-defects in GICs, such as crazing and voids, are responsible for the weak mechanical properties.

Turning our attention towards the aqueous solutions used in this study, distilled water and artificial saliva, higher flexural strength was measured in both solutions under coating conditions. Within a coating condition, the storage media were proved to have only scattered significant effect on the physical properties measured. Our finding of a low effect of the storage agent on the measured mechanical properties is in agreement with Nicholson and Wilson, who stated that there are no statistically significant differences of strengths measured between the storage media, among others deionized water and artificial saliva, concerning the time period from 24 h to 30 days [30]. Focusing on the condition of storage time, the different materials evolved distinguished behaviour of mechanical strength. In view of the measured micro-mechanical properties, VH and *E* [Table 4 (a and b)], these properties are likely to depend stronger on the composition of the selected materials than the measured macro-mechanical properties. Many studies [22, 31-34] discussed this assumption as they described the influence of the chemical composition, concentration and molecular weight of the polycarboxylic acid, the glass structure and the power/liquid ratio.

In contrast to the flexural strength, the effect of resin coating on the micro-mechanical properties was mostly not significant. A particular behaviour was observed for Chem-Fil Rock. Though reaching significant higher macromechanical properties when compared to the other GICs, the micro-mechanical properties were the lowest. The reason for lower micro-mechanical properties must be searched in the filler size and morphology [23] of the glass particles in ChemFil Rock. Analysing the micro-mechanical attributes of GIC and the effect of coating, it has to be ascertained that surface protection in the form of resin coating did not show the expected effect of a significant improvement in hardness and indentation modulus.

Conclusion

The new GIC might represent a promising approach of granting GICs' higher longevity, which turns it into a more favourable filling material in class I and II cavities. In contrast to the other materials tested, it showed no crack propagation and visible surface irregularities and thus improved macromechanical characteristics, but lower micro-mechanical properties. The micro-mechanical properties of the analysed GICs did neither benefit from the new zinc formulation nor from resin coating. A resin coating is nevertheless a valuable support for GIC fillings, since it led to significant improvements in flexural strength. This statement is valid also for ChemFil Rock, contrary to the manufacture indication. The impact of storage agent and duration on the measured properties was low.

Conflict of interest The authors declare that they have no conflicts of interest.

References

- Glasspoole EA, Erickson RL, Davidson CL (2002) Effect of surface treatments on the bond strength of glass ionomers to enamel. Dent Mater 18:454–462
- 2. Hunt PR (1994) Glass ionomers: the next generation. A summary of the current situation. J Esthet Dent 6:192–194
- Mount GJ (1994) Buonocore Memorial Lecture. Glass-ionomer cements: past, present and future. Oper Dent 19:82–90
- 4. Wiegand A, Buchalla W, Attin T (2007) Review on fluoridereleasing restorative materials—fluoride release and uptake

characteristics, antibacterial activity and influence on caries formation. Dent Mater 23:343-362

- Naasan MA, Watson TF (1998) Conventional glass ionomers as posterior restorations. A status report for the American Journal of Dentistry. Am J Dent 11:36–45
- Ilie N, Hickel R, Valceanu A, Huth K (2012) Fracture toughness of dental restorative materials. Clin Oral Invest 16:489–498
- Crisp S, Wilson AD (1974) Reactions in glass ionomer cements: I. Decomposition of the powder. J Dent Res 53:1408–1413
- Earl MS, Mount GJ, Hume WR (1989) The effect of varnishes and other surface treatments on water movement across the glass ionomer cement surface. II. Aust Dent J 34:326–329
- Hotta M, Hirukawa H, Yamamoto K (1992) Effect of coating materials on restorative glass-ionomer cement surface. Oper Dent 17:57–61
- Ribeiro AP, Serra MC, Paulillo LA, Rodrigues Junior AL (1999) Effectiveness of surface protection for resin-modified glassionomer materials. Quintessence Int 30:427–431
- Watson T, Banerjee A (1993) Effectiveness of glass-ionomer surface protection treatments: a scanning optical microscope study. Eur J Prosthodont Restor Dent 2:85–90
- DENTSPLY Wissenschaftliches Kompendium ChemFil Rock. www.dentsply.de. Accessed 14 Feb 2011
- 13. Wilson ADMJ (1988) Glass ionomer cement. Quintessence Publishers, Chicago
- Crisp S, Lewis BG, Wilson AD (1976) Characterization of glassionomer cements 1. Long term hardness and compressive strength. J Dent 4:162–166
- Mitra SB, Kedrowski BL (1994) Long-term mechanical properties of glass ionomers. Dent Mater 10:78–82
- Pearson GJ, Atkinson AS (1991) Long-term flexural strength, of glass ionomer cements. Biomaterials 12:658–660
- Crisp S, Lewis BG, Wilson AD (1979) Characterization of glassionomer cements 5. The effect of the tartaric acid concentration in the liquid component. J Dent 7:304–312
- Wasson EA, Nicholson JW (1991) Studies on the setting chemistry of glass-ionomer cements. Clin Mater 7:289–293
- Mount GJ, Hume WR (1998) Preservation and restoration of tooth structure. Harcourt Brace & Company Ltd/Mosby International, London
- Gemalmaz D, Yoruc B, Ozcan M, Alkumru HN (1998) Effect of early water contact on solubility of glass ionomer luting cements. J Prosthet Dent 80:474–478

- Cattani-Lorente MA, Godin C, Meyer JM (1993) Early strength of glass ionomer cements. Dent Mater 9:57–62
- 22. Prosser HJ, Powis DR, Wilson AD (1986) Glass-ionomer cements of improved flexural strength. J Dent Res 65:146–148
- Xie D, Brantley WA, Culbertson BM, Wang G (2000) Mechanical properties and microstructures of glass-ionomer cements. Dent Mater 16:129–138
- De Barra E, Hill RG (2000) Influence of glass composition on the properties of glass polyalkenoate cements. Part III: influence of fluorite content. Biomaterials 21:563–569
- Darling M, Hill R (1994) Novel polyalkenoate (glass-ionomer) dental cements based on zinc silicate glasses. Biomaterials 15:299–306
- Boyd D, Towler M, Law R, Hill R (2006) An investigation into the structure and reactivity of calcium-zinc-silicate ionomer glasses using MAS-NMR spectroscopy. J Mater Sci Mater Med 17:397– 402
- Boyd D, Towler M (2005) The processing, mechanical properties and bioactivity of zinc based glass ionomer cements. J Mater Sci Mater Med 16:843–850
- Moshaverinia A, Roohpour N, Darr JA, Rehman IU (2009) Synthesis and characterization of a novel *N*-vinylcaprolactam-containing acrylic acid terpolymer for applications in glass-ionomer dental cements. Acta Biomater 5:2101–2108
- Prentice LH, Tyas MJ, Burrow MF (2005) The effect of particle size distribution on an experimental glass-ionomer cement. Dent Mater 21:505–510
- Nicholson JW, Wilson AD (2000) The effect of storage in aqueous solutions on glass-ionomer and zinc polycarboxylate dental cements. J Mater Sci Mater Med 11:357–360
- Wilson AD, Hill RG, Warrens CP, Lewis BG (1989) The influence of polyacid molecular weight on some properties of glass-ionomer cements. J Dent Res 68:89–94
- Crisp S, Lewis BG, Wilson AD (1976) Characterization of glassionomer cements: 2. Effect of the powder: liquid ratio on the physical properties. J Dent 4:287–290
- Crisp S, Lewis BG, Wilson AD (1977) Characterization of glassionomer cements. 3. Effect of polyacid concentration on the physical properties. J Dent 5:51–56
- Griffin SG, Hill RG (1999) Influence of glass composition on the properties of glass polyalkenoate cements. Part I: influence of aluminium to silicon ratio. Biomaterials 20:1579–1586

Copyright of Clinical Oral Investigations is the property of Springer Science & Business Media B.V. and its content may not be copied or emailed to multiple sites or posted to a listserv without the copyright holder's express written permission. However, users may print, download, or email articles for individual use.